

TRANSFORMATION, COLONY PICK UP, MINIPREP, DIGESTION

Transformation of Bacteria by heatshock method

Note, it is not correct to say “transformation of plasmid”
TA will do up to 2 for you

1. Prepare ice in ice bucket
2. Thaw competent cell (bacteria) on ice. Always keep on ice.
3. Transfer 90 μ l bacteria in precooled 15 ml falcon tubes.
4. Add ~10 μ l (whatever left in tube) of ligation product
5. Incubate on ice for 30 min.
6. Warm at 42 °C for 45 sec and return on ice.
7. Add 1 ml of LB.
8. Incubate at 37 °C for 1 hour with shaking.
9. While waiting, leave the kanamycin plates (one for each reaction) at 37 °c with the lid half open.
10. Apply 200 μ l on to kanamycin plate.
11. Add 10-20 glass beads
12. Shake the plate horizontally and let the beads roll around to spread the bacteria.
13. Incubate at 37 °C for overnight. If the plate is wet, leave the lid open until it dries.

Pick-up colonies

1. Prepare terrific-broth (TB medium)+50 μ g/ml kanamycin in 15 ml culture tubes (approx. 1.5 ml/tube).
2. Identify a well-separated colony.
3. Touch with pipette tip or sterile forceps and drop it in the tube.
4. Shake in 37 °C incubator for O/N.

Miniprep

1. Grow minipreps in TB medium+50 μ g/ml kanamycin, from single colonies in 15 ml culture tubes (approx. 1.5 ml/tube). Each 4 per group.
2. Simply pour bacteria from the Falcon tubes into the epp tubes. Pour 1-1.5 ml. Don't spill b/c it will cross contaminate. Keep original culture tubes.
3. Spin epp. tubes at max. speed for 15 sec.
 - Place the epp tubes in the centrifuge with the cap opening facing down.
 - Speed maximum is at 13.0
4. Open the tubes carefully, watching for cross contamination. Dump supernatant by pouring it into the waste beaker from a sort of high height (you don't want any contamination). Make sure not to mix the contents of different tubes.
5. add 200 μ l of solution 1 to bacterial pellet.
 - Solution 1 contains the following:

4.5 grs glucose
12.5 ml of Tris pH 8
10 ml EDTA pH 8
RNase

fill with distilled water up to 500 ml.

- EDTA eliminates Ca and Mg so DNase won't be active and therefore will not degrade DNA. Thus, the DNA is stable.
- Tris is a buffer.

6. Close eppendorf tubes and vortex them (until pellet is dissolved).
7. Add 200 μ l of solution 2 and close tubes. Invert each tube gently to mix the contents well.
 - Solution 2 contains the following:
 - 10 ml of NaOH 10N
 - 50 ml of 10% SDS
 - fill up to 500 ml with water
 - SDS is a detergent that lyses cells
 - NaOH is a strong base that works with the strong acid in Solution 3 to cause genomic DNA to precipitate leaving the plasmid DNA in solution.
8. Open tubes and add 200 μ l of solution 3. Mix vigorously by shaking with your hand.
 - Solution 3 contains the following:
 - 89 grs of potassium acetate
 - 58 ml of acetic acid
 - fill up to 500 ml with distilled water
9. Spin at max. speed for 3-5 minutes at RT.
 - The plasmid DNA remains in the supernatant. The genomic DNA is the pellet.
10. During the spin, prepare new epp. tubes, and add 1 ml of isopropanol to each tube.
 - DNA precipitates with isopropanol. The minimum volume necessary for this is 0.7x your solution.
11. Take 500 μ l of the supernatant from the spin tubes and add to the tubes containing 1 ml of isopropanol.
 - Try to get just the liquid, and leave behind the gunk from the tube. Mix vigorously by hand.
 - Pipet gingerly to avoid the pellet. Tilt the tube nearly horizontally and run the pipet down the wall of the tube.
 - The DNA is very stable in alcohol.
12. Spin max. speed RT for 5-10 minutes
13. Dump supernatant. (See notes for dumping the supernatant for step 4).
14. Add ~500 μ l of 70% ethanol. Spin for 1 minute.
15. Dump ethanol. Keep tubes inverted over paper towel to drain the remnants of ethanol.

-Make sure to set the tubes far enough apart to avoid cross contamination.

16. After DNA is dry, resuspend in 50 μ l 10 mM Tris HCl (pH 8.0).

Digestion of Miniprep DNA

1. Prepare digestion mix:

Reagent	Quantity 5X (μ L)	Quantity 1X (μ L)
Buffer 2	5	1
EcoRI	2.5	0.5
XhoI	2.5	0.5
(bovine serum albumin) BSA	5	1
dH ₂ O	20	4

Measure out each reagent to the 5X amount. Then use this solution to aliquot one-fifth of the solution (1X) into each tube for your digestion.

NOTE: The enzymes are not stable at room temperature: keep them out for as short as it is possible, or put them in a cooler, or on ice.

-Centrifuge the mix about 10 seconds in its epp tube labeled “digestion mix”

2. Prepare tubes for digestion: aliquot 7 μ l of the digestion mix/tube

3. Vortex and add 3 μ l of the DNA into each 7 μ l digestion tube. Mix well by pipetting.

4. Digest the 4 DNA's that you purified.

5. Incubate at 37C for at least 40 minutes (some enzymes require longer).

6. Run entire sample for each tube on 0.8% agarose gel along with PCR product and vector DNA you prepared last week and along with the DNA ladder.

-Do not forget to add BPB to samples if you have not done so previously.